See discussions, stats, and author profiles for this publication at: https://www.researchgate.net/publication/309445272

Drug-eluting coating of ginsenoside Rg1 and Re incorporated poly(lactic-co-glycolic acid) on stainless steel 316L...

Article in International Journal of Pharmaceutics · October 2016

DOI: 10.1016/j.ijpharm.2016.10.056

FIONS		reads 32	
thors	, including:		
	Mun Fai Loke		Syafiqah Saidin
dir.	National University of Singapore		Universiti Teknologi Malaysia
	100 PUBLICATIONS 218 CITATIONS		14 PUBLICATIONS 81 CITATIONS
	SEE PROFILE		SEE PROFILE
	Hendra Hermawan		
	Laval University		
	84 PUBLICATIONS 1,046 CITATIONS		
	SEE PROFILE		

Some of the authors of this publication are also working on these related projects:




Proteomics and Genomics approaches for studying obesity related gastric cancer: Facilitate the development of new therapy for Immunocompetence View project



High strength biodegradable metals for temporary biomedical implants View project

All content following this page was uploaded by Hendra Hermawan on 29 October 2016.

The user has requested enhancement of the downloaded file. All in-text references <u>underlined in blue</u> are added to the original document and are linked to publications on ResearchGate, letting you access and read them immediately.



http://dx.doi.org/10.1016/j.ijpharm.2016.10.056

Get rights and content

Abstract

Active ingredients of ginsenoside, Rg1 and Re, are able to inhibit the proliferation of vascular smooth muscle cells and promote the growth of vascular endothelial cells. These capabilities are of interest for developing a novel drug-eluting stent to potentially solve the current problem of late-stent thrombosis and poor endotheliazation. Therefore, this study was aimed to incorporate ginsenoside into degradable coating of poly(lacticco-glycolic acid) (PLGA). Drug mixture composed of ginseng extract and 10% to 50% of PLGA (xPLGA/g) was coated on electropolished stainless steel 316L substrate by using a dip coating technique. The coating was characterized principally by using attenuated total reflectance-Fourier transform infrared spectroscopy, scanning electron microscopy and contact angle analysis, while the drug release profile of ginsenosides Rg1 and Re was determined by using mass spectrometry at a one month immersion period. Full and homogenous coating coverage with acceptable wettability was found on the 30PLGA/g specimen. All specimens underwent initial burst release dependent on their composition. The 30PLGA/g and 50PLGA/g specimens demonstrated a controlled drug release profile having a combination of diffusion- and swelling-controlled mechanisms of PLGA. The study suggests that the 30PLGA/g coated specimen expresses an optimum composition